

Geon™ Vinyl Flexible HC38000G Flexible Polyvinyl Chloride

Key Characteristics

Product Description

Geon HC 3180G is a non-phthalate, gamma-stable flexible vinyl compound made from FDA sanctioned ingredients. HC 3180G can be used for extrusion and injection molding healthcare applications.

This product, with color designation as "Clear 9883",has obtained Certification under Class VI U.S. Pharmacopea for systemic toxicity and intractaneus studies.

General			
Material Status	Commercial: Active		
Regional Availability	 Africa & Middle East Asia Pacific	EuropeLatin America	North America
Features	BPA FreeEthylene Oxide SterilizableGood Processability	High ClarityNon-Phthalate PlasticizerRadiation (Gamma) Resistant	Radiation Sterilizable
Uses	ConnectorsDisposable Hospital GoodsHospital Goods	 Medical Devices Medical/Healthcare Applications Personal Care 	• Tubing
Agency Ratings	 FDA Unspecified Rating 		
Appearance	 Clear - Blue Tint 		
Forms	• Pellets		
Processing Method	 Extrusion 	Injection Molding	

Technical Properties 1

Physical	Typical Value (English)	Typical Value (SI)	Test Method
Density / Specific Gravity	1.23	1.23	ASTM D792
Mechanical	Typical Value (English)	Typical Value (SI)	Test Method
Tensile Strength ² (100% Strain)	1140 psi	7.86 MPa	ASTM D638
Tensile Strength ² (Break)	2700 psi	18.6 MPa	ASTM D638
Tensile Elongation ² (Break)	440 %	440 %	ASTM D638
Hardness	Typical Value (English)	Typical Value (SI)	Test Method
Durometer Hardness (Shore A, 15 sec)	80	80	ASTM D2240
Thermal	Typical Value (English)	Typical Value (SI)	Test Method
Brittleness Temperature	-13.0°F	-25.0 °C	ASTM D746
Additional Information			

Additional Information

NOTE: All ingredients used in the formulation of this compound are listed in the United States FDA Code of Federal Regulations, Title 21. The only restriction is that the plasticizer is "For foods of high water content only". This statement is only an overview of the FDA listings for this product. It is the customer's responsibility to ensure that the compound is suitable for its intended application.

This product is not recommended for long-term storage of blood or blood components. Please contact PolyOne for more information.

Notes

¹ Typical values are not to be construed as specifications.

Rev: 2016-01-26 Page: 1 of 2

² Type IV, 20 in/min (510 mm/min)